

1.1 PRODUCT DESCRIPTION OVERVIEW

The Abbott Aeroset® System is an automated, random-access clinical chemistry system. The instrument uses a dual-pipetting system with a maximum throughput of up to 2000 tests/hour. It can be configured with up to 100 different assay files utilizing photometric and potentiometric methodologies. The potentiometric assays (Na^+ , K^+ , and Cl^-) are measured using an Integrated Chip Technology™ (ICT™). The Abbott Aeroset System will process routine and STAT samples using independent Carousel and FastTrack™ Samplers.

1.2 HOW TO USE THIS MANUAL

The purpose of this manual is to provide information useful for servicing the Abbott Aeroset® System and is composed of seven (7) chapters.

Chapter 1. General Data

Contains a product overview, information on manual usage, manual revision marks, accident prevention symbols, and system specifications.

Chapter 2. Troubleshooting

Contains procedures to be used by Field Service Engineers/Representatives in identifying and isolating problems. The procedures are presented in probable cause and resolution format.

Chapter 3. Parts Lists

Provides the Supplemental Tools and Supplies list, exploded view drawings of all field spared parts (and their configuration) which are indexed to the item number, and description on the spare Parts Lists (PL).

Chapter 4. Removal and Replacement

Contains Removal & Replacement procedures that are indexed by number to related Parts Lists in Chapter 3. For example, RR-A1.5 is the Removal & Replacement procedure for Parts List PL-A1, item #5.

Chapter 5. Verification Procedures

Contains adjustments, calibrations, checks, tests, and additional procedures required to verify instrument operation after repairs are complete. Verification Procedures (VPs) are also used to assist in troubleshooting.

Chapter 6. PM/Total Call

Contains preventive maintenance (PM) procedures, Total Call procedure, and PM and Total Call checklists for the instrument.

Chapter 7. Installation

Outlines the responsibilities of the Field Service Engineer/ Representative (FSE/FSR) in setup of the Abbott Aeroset System from pre-installation to installation and relocation.

1.2 HOW TO USE THIS MANUAL, CONT.

Manual Revision Marks

Abbott Aeroset® Service Manual revision pages keep the manual up to date with configuration changes and servicing techniques. The actual changes will be identified as follows:

Text Revisions

- Black bar in left margin for significant changes that impact servicing the instrument.

Art Revisions

Black bar left of art number for significant changes that impact servicing the instrument.
(Figure 5 - 1. Revision Marks)

Page Revisions

If anything on a page changes, the text part number will be changed by incrementing the dash (-) number for that page (e.g., from -101 to -102).

A new title page with the revision log will be sent with each change package. This list will contain the page number of each changed or added page along with the revision dash (-) number of that page. Pages not listed are original pages.




Figure 1-1. Revision Marks

1.2 HOW TO USE THIS MANUAL, CONT.




Identification of Unique Items

If information, a procedure or a spare part is unique to United States Field Service only, then **(US)** is noted beside the item. If the item is unique to International Field Service only, **(INTL)** is noted beside item. If information, a procedure, or spare part is unique to another country, the country is listed in parenthesis beside that item.

Example:

<u>ITEM</u>	<u>PART NO.</u>	<u>DESCRIPTION</u>
11	2E12-01	Keyboard (English)
	2E08-01	Keyboard (French)
	2E09-01	Keyboard (Spanish)

Special Symbols



SYMBOL	DEFINITION
GND	Chassis GND (Ground) unless otherwise stated.
	Removal & Replacement Procedure (RR).
	Verification Procedure (VP).
	<p>ELECTROSTATIC SENSITIVITY</p> <p>The Electrostatic Discharge symbol identifies an activity or area where the operator must wear a ground strap while servicing the system.</p> <p>Components sensitive to electrostatic discharge such as printed circuit boards, computer components, etc. may be damaged.</p>

1.2 HOW TO USE THIS MANUAL, CONT.

TSB/ISA Symbols


Two symbols are used to show areas or sections in the service manual which have been affected by a TSB or an ISA. Refer to actual TSB or ISA for additional information.

These symbols are used to show a particular part or area which has been, or has not been, modified by the TSB number specified within the symbol.



TSB 125-25 Installed TSB 125-25 Not Installed

This symbol is used to identify the ISA containing additional information about the part or area.



1.3 HAZARD TYPES

Introduction

Operation, maintenance, and servicing of the Abbott Aeroset® Systems may expose individuals to potential safety and health hazards. All work must be performed according to procedures described in the Abbott Aeroset Operations or Service Manuals or as directed by an Abbott Technical Support Representative. This section describes the types and locations of potential hazards that could cause physical harm to service personnel.

Hazard Signal Word Definitions and Icons

Hazard Signal Words

DANGER

Denotes an Immediate hazard which, if not avoided, could result in serious injury or death. This signal word represents the highest level of any hazardous situation.

WARNING

Denotes a hazard which could result in moderate or serious personal injury.

CAUTION





Denotes potential hazards that could result in minor injury. Also used for conditions or activities that could interfere with the proper functioning or performance of the instrument.

NOTE: *Denotes general operator or service information that will not impact performance or results if it is not observed.*

1.3 HAZARD TYPES, CONT.




Hazard Icons

Safety hazard icons are used in this manual and on the instrument labels to identify potentially dangerous conditions or situations. In this manual and on some instrument labels, text accompanies the safety icon to describe the hazard. For other instrument labels, the operator or service technician is to refer to the manual for specific information. All users of the instrument must be familiar with the following messages:

ICON	RELATED TEXT
	DESCRIPTION
	<p>WARNING: POTENTIAL BIOHAZARD.</p> <p>Identifies an activity or area where the operator may be exposed to potentially infectious materials or substances.</p>
	<p>DANGER: HIGH VOLTAGE.</p> <p>Identifies high voltage areas over 600 volts.</p>
	<p>WARNING: ELECTRICAL SHOCK HAZARD.</p> <p>Identifies the possibility of electrical shock if procedural or engineering controls are not observed.</p>
	<p>CAUTION: HOT SURFACE</p> <p>Identifies a part or component that requires cooling before touching.</p>

1.3 HAZARD TYPES, CONT.

Hazard Icons, Cont.

ICON	RELATED TEXT
	DESCRIPTION
	<p>CAUTION: CLASS 2 LASER PRODUCT</p> <p>Warns against direct viewing into the Bar Code Reader laser beam or reflections of the beam from reflective surfaces.</p>
	<p>CAUTION: (ASSOCIATED TEXT)</p> <p>Identifies an activity or an area that may present a physical, mechanical, or chemical hazard. The user must be aware, alert, and cautious to prevent physical injury.</p>
	<p>WARNING: (ASSOCIATED TEXT)</p> <p>See above. A more severe degree of hazard than above.</p>

1.4 BIOHAZARD SAFETY



WARNING: POTENTIAL BIOHAZARD

Consider all clinical specimens, reagent, controls, and calibrators, etc. that contain human blood and surfaces or components that have come into contact with human blood as potentially infectious. Wear gloves, lab coats, and safety glasses, and follow other biosafety practices as specified in OSHA Bloodborne Pathogen Rule (29 CFR Part 1910.1030) or other equivalent biosafety procedures.

Handling Spills

Clean up spills of potentially infectious materials according to established biosafety practices. Use the following generally accepted procedure for cleaning such spills.

1. Absorb the spilled material with absorbent towels or other absorbent materials.
2. Wipe the area with detergent solution

3. Wipe the area with an appropriate disinfectant such as 2% chlorine bleach solution or a 0.1% sodium hypochlorite solution.

Instrument Decontamination

Instruments must be decontaminated prior to servicing or shipment. Complete the following procedure unless an Abbott Technical Support Representative provides other instructions

1. Cycle 10% chlorine bleach solution through the fluid pathway (including the probe) that contacted blood or products containing blood. Allow the bleach solution to remain in the system for a minimum of ten minutes.
2. Cycle water or buffer solution through the system to rinse the affected components.
3. Remove any liquid material from the instrument.
4. Rinse waste containers with the chlorine bleach solution or other suitable disinfectants before shipping or disposal.

1.4 BIOHAZARD SAFETY, CONT.

5. Wipe the surface of each instrument with a detergent solution followed by an appropriate disinfectant such as 0.1% sodium hypochlorite (2% chlorine bleach solution). Allow to air dry for a minimum of ten minutes.

Handling Waste

Dispose of all clinical specimens, reagents, controls, calibrators, and other disposables that may be contaminated according to local, state, federal, and country regulations governing the treatment of regulated medical waste.

Sharps

Sharps, such as contaminated probes, should be placed in appropriately marked, puncture-resistant containers before treatment and disposal.

1.5 LASER SAFETY

CAUTION: CLASS 2 LASER PRODUCT

Denotes low-power, visible-light lasers or laser systems. This class of laser do not normally present a hazard to the eye because of normal human aversion response such as blinking, eye-movement, and the like. However, like many conventional light sources, they may present some potential for hazard if viewed directly for extended periods of time.

During normal instrument operation, the inner protective covers (and access door) are to remain in place to prevent unintentional laser light exposure. The inner protective covers (and access doors) may be removed by qualified personnel during servicing.

Do not remove or alter the inner protective cover laser warning labels. They must remain legible.

The Class 2 Laser label on the instrument is shown below in Figure 5 - 2. Laser Caution Label. The label consists of black lettering on a yellow background.



Figure 5-2. Laser Caution Label

1.6 MECHANICAL/PHYSICAL SAFETY



WARNING: PUNCTURE HAZARD

The probe is sharp. Avoid contact with the tip of the probe. Use protective gloves if necessary.



CAUTION: HOT SURFACE

Specific components and parts are very hot. Allow the part or component to cool sufficiently (generally 10 minutes) before servicing or replacement. Use temperature resistant gloves if required.



WARNING: LIFTING HAZARD

The Abbott Aeroset® System weighs 1500 pounds (680 kg). Obtain assistance when moving and/or use appropriate lifting devices.

1.7 ELECTRICAL SAFETY



WARNING: ELECTRICAL SHOCK HAZARD

Turn off the power to the instrument and disconnect the power cord before replacing fuses, printed circuit boards, and other electrical parts and components. Replace only the fuses that are externally accessible and labeled. Use only replacement fuses of the specified type and electrical rating.

Keep surfaces dry around the instrument and electrical connectors.

Use approved power cords and electrical accessories only, such as those supplied with the instrument to protect against electrical shock. Connect power cords to properly grounded outlets only.

Do not disconnect any electrical connections while the power is ON. Follow instructions for correctly powering-down the instrument and all connected equipment before performing service or maintenance.

Shutdown Procedure

Prior to turning system power OFF, perform the following procedure.

1. From the DATABASE screen:
 - SHUTDOWN
2. SHUTDOWN CONFIGURATION screen is displayed.
3. Check boxes for procedures to be performed.
 - a. For 24-hour operation, select:
 - Clear Database, All Records
 - b. For partial-day operation, select:
 - Change Water in Bath
 - Wash Cuvettes
 - Wash Probes with - 1% Alkaline Solution
 - Wash ICT™ with ICT™ Cleaning Fluid and ICT™ Reference Solution
 - Drain and Fill IRef Cup
 - Clear Database, All Records
4. Check Acid and Alkaline volumes and replace if necessary.

1.7 ELECTRICAL SAFETY, CONT.

5. When SHUTDOWN CONFIGURATION is complete, select:

■ SHUTDOWN

If the Rotary Power Control Switch is set to AUTO, system power will be shutdown after the shutdown process.

High Voltage Hazard



DANGER: HIGH VOLTAGE

Ensure the power to the instrument is turned OFF. A high voltage charge may remain on electrical components such as laser connectors or power supplies with the system power off. Use an electrically insulated tool to disconnect the charged part or component and short the charge (both male pins) to the instrument chassis.

1.8 CHEMICAL SAFETY



CAUTION: CHEMICAL HAZARD

Wear safety glasses, chemical resistant gloves, and a lab coat when handling the Abbott solutions. These solutions may be potentially harmful. Refer to the Material Safety Data Sheet (MSDS) or the package insert for specific safety information. In case of contact with the skin or eyes, flush with water for at least 15 minutes. If irritation persists or signs of toxicity occur from exposure, seek medical attention immediately.

1.9 SYSTEM SPECIFICATIONS

This section contains instrument dimensions, space, computer and interface specifications, electrical, environmental, printer, and optical specifications and capacities requirements.

Physical Dimensions

Analyzer

Width:	74.5" (189 cm)
Depth:	44.1" (112 cm)
Height:	43.5" (110.5 cm)
Weight:	1500 lbs. (680 kg)

System Control Center (SCC)

Width:	24" (61 cm)
Depth:	18" (46 cm)
Height:	42.5" (108 cm)
(with monitor)	60.5" (154 cm)
Weight:	190 lbs. (86.5 kg) approx.

Clearances

Right:	31.5" (80 cm)
Left:	31.5" (80 cm)
Rear:	20" (51 cm)
Above:	20" (51 cm)
Front:	36" (92 cm)

For more detailed system dimension information, see [Chapter 7, Installation](#).

1.9 SYSTEM SPECIFICATIONS, CONT.

Environmental Requirements

NOTE: *Instrument For Indoor Use Only*

Location: Flat, level surface. No direct sunlight or drafts. Remove from sources of direct heat and moisture. Do Not place next to a heat or vibration-generating device.

Temperature -
Operational: Minimum 15° C (59° F)
Maximum 32° C (89° F)

Storage
Temperature: Minimum 5° C (41° F)
Maximum 50° C (122° F)

Humidity: 15 - 80%
(non-condensing)

BTU Output: 6551 BTU/Hr.

Operating
Altitude: 6600 ft. (2000 m)

Computer and Interface Specifications

Processor: Pentium®, 133 MHz, 16 MB RAM

Hard Drive: 2.1 GB, EIDE bus interface to CPU

Disk Drive: 3.5", 1.44 Mb, IDE bus interface to CPU

Operator
Interface Display: Full color VGA CRT, 17" diagonal (43cm)

Keyboard: 101 Key, custom, IBM® compatible

Printer: 80-column, parallel-port, Centronix® connection

Host Interface
(System): Bi-directional, RS-232 serial communication port, 1200, 2400, 4800, or 9600 baud

1.9 SYSTEM SPECIFICATIONS, CONT.

Electrical Specifications

Required Measurement at 220V Outlet

U.S.	BREAKER CLOSED	BREAKER OPEN
Line 1 to Line 2 =	200 - 240 VAC	≤ 0.5 VAC
Line 1 to GND =	90 - 132 VAC	≤ 0.5 VAC
Line 2 to GND =	90 - 132 VAC	≤ 0.5 VAC
GND to Conduit =	< 0.5 VAC	≤ 0.5 VAC

INTERNATIONAL	BREAKER CLOSED	BREAKER OPEN
Line to Neutral =	180 - 264 VAC	≤ 0.5 VAC
Line to Earth =	180 - 264 VAC	≤ 0.5 VAC
Neutral to Earth =	0.5 VAC	≤ 0.5 VAC
Neutral to Case =	< 0.5 VAC	≤ 0.5 VAC

Analyzer

- Single phase, 180 - 264 VAC.
- 50/60 Hz ± 1%.
- 20 Amps - Nominal Operation.

Dedicated Power Line

- Separate Circuit breaker/Fuse.
- No other electrical devices can be supplied from dedicated circuit breaker and power lines.

Conduit or BMx (Flexible Metal Conduit)

- No other electrical wires are present

Receptacle (U.S./Canada)

- The receptacle shall be a NEMA L6-30R

1.9 SYSTEM SPECIFICATIONS, CONT.

Electrical Specifications, cont.

Power Connection

Receptacle Mounting	The electrical box containing the receptacle shall be mounted within 12-ft. (3.7m) of the left side of the instrument.
Power Cord (U.S./Canada)	UL/CSA-approved, SJT type, 12 AWG, 250VAC.
Power Cord (Outside U.S./ Canada)	IEC-309 approved for at least 30A, 250VAC and be of a type meeting national electrical code of the country of final destination.
Analyzer Connector: (U.S./Canada)	The service connector shall be a NEMA L6-30P.
Analyzer connector: (outside U.S./ Canada)	IEC-309 approved for at least 30A, 250VAC and be of a type meeting national electrical code of country of final destination.

Printer

Power Requirements

110VAC areas:	108 - 132VAC
220VAC areas:	198 - 264VAC
Frequency:	50/60 Hz \pm 1%
Service Conn. (U.S./Canada)	IEC 320, UL/CSA-rated
Service Connector: (Outside U.S./ Canada)	Plug on power cord, mating with utility service must be rated and IEC-approved for at least 120VAC and 0.4A (48 VA) and be of a type meeting national electrical code of country of final destination. The plug is to be attached at final destination if other than supplied plug.

System Control Center

Power Requirements:
System power supplied from analyzer.

1.9 SYSTEM SPECIFICATIONS, CONT.**Sample Bar Code Specifications**

BAR CODE READER	SUPPORTED SYMBOLOGY
Sample	Codabar Code 3 of 9 Code I 2 of 5 Code 128, Subset A, B, and C 4 - 20 Characters, No Checksum or Start/Stop Characters
Reagent	2D Labels PDF - 417 1D Labels Code I 2 of 5 Code 128 Codabar Code 3 of 9

Measurement MethodsType

Photometric

Potentiometric

Throughput

Photometric and ICT™: Up to 2000/hr.

Photometric only: Up to 1600/hr.

ICT only: Up to 600/hr

1.9 SYSTEM SPECIFICATIONS, CONT.**Optical Specifications**

Lamp Type:	Tungsten - Halogen
Wavelength Range:	340 - 804 nm
Wavelengths Measured: (in nanometers)	340, 380, 404, 412, 444, 476, 500, 524, 548, 572, 604, 628, 660, 700, 748, and 804 nm
Photometric Range:	0.1 to 3.0 Abs (converted to 10 mm light path length)
Light Path Length:	5 mm
Linearity:	Within $\pm 2\%$ @ 2.0 Abs

CapacitiesData Storage

Data Base Records:	8000
QC:	150/Level
Error Log:	200

Carousels

Reaction:	330 Reaction Cuvettes
Cal/Control:	45 Positions
Sample:	30 Positions
Stat:	1 Position
Reagent:	56 Positions/Carousel
Seg. A:	12 Positions
Seg. B:	12 Positions
Seg. C:	12 Positions
Seg. D:	20 Positions

1.9 SYSTEM SPECIFICATIONS, CONT.**Sample Carrier & Sample Carrier Tray**

Sample Carrier:	5 Samples per Carrier
Sample Carrier Tray:	10 Sample Carriers per tray
FastTrack™ Capacity:	200 samples (4 Sample Carriers)

Sample Vessels

<u>Sample Vessel</u>	<u>Ext. Diam (mm) Hgth (mm)</u>
10-mL (16 mm) tube:	15.4 ± 0.3100
7-mL (16 mm) tube:	15.4 ± 0.375
7-mL (13 mm) tube:	12.4 ± 0.2100
5-mL (13 mm) tube:	12.4 ± 0.175
Sample Cup (2 mm)	

Reagent Containers

Large:	90 mL
Small:	50 mL (requires adapter)
100 mL:	100 mL
20 mL:	Round (requires adapter)

Bulk SolutionsBulk Solution

ICT™ Ref. Sol.:	1 L
Alkaline Wash:	500 mL
Acid Wash:	500 mL

Method Volumes

Sample Dispense:	2-35 µL/Test
Sample Aspiration:	Sum of sample Dispense
Sample ICT™	
Dispense:	15 µL
Reagent Dispense:	20-345 µL/test
Reagent Aspiration:	Sum of Reagent Dispense

Liquid Waste Container: Approx. 8 L

1.9 SYSTEM SPECIFICATIONS, CONT.**Fluidics**

TUBING TYPE	ID	OD	LGTH.
Deionized Water High Pressure Hose	9mm	5mm	≤5 M
Low Concentrate Waste High Pressure Hose	15mm	23mm	≤5 M
Overflow High Pressure Hose	12mm	18mm	≤5 M
Bath Water Drainage High Pressure Hose	12mm	18mm	≤5 M
High Concentrate Waste Silicon Rubber	6mm	12mm	≤5 M

Drainage

Water Drainage (Total During Analysis)	47.5 L/hr
Maximum Drainage	70 L/hr
High Concentrate Drainage	2.5 L/hr